

pH Buffers and Conductivity Standards

For calibrating pH and Conductivity meters

Features:

pH Buffers

- Easy to use, already made fresh pH buffers for calibration
- Helpful table shows what the pH value should be vs. the temperature of the buffer solution measured
- Rinsing solution used for cleaning off the pH electrode after taking the next measurement from a different solution to prevent contamination
- Convenient pouch eliminates the need for a cup to store buffer or rinsing solution
- PH103 contains 6 pouches of each buffer (pH 4.00, 7.00, and 10.00), plus 2 Rinsing solutions
- Pint size pH buffer solutions available (PH4-P, PH7-P, and PH10-P)
- 650470 contains 3 vials (10 capsules each) of 4, 7, and 10 pH dry powder buffers and one dispenser vial of color key preservative. Mix with 100mL DI water.

Conductivity Standards

- Convenient pint size conductivity standards calibrate your conductivity meter to optimize its accuracy
- Available standards: 84 μ S, 1413 μ S, and 12880 μ S



PH103 includes 6 pouches each of pH 4, pH 7, and pH 10 Buffers plus 2 Rinse solutions



650470 pH 4, 7, and 10 buffer capsules with color key preservative



Pint size pH buffers for pH 4, 7, and 10 available (PH4-P, PH7-P, PH10-P)



EC-84-P, EC-1413-P, and EC-12880-P pint size conductivity standards

Ordering Information:

- PH103** Tripak Buffer Pouches (6 ea 4,7,10 pH & 2 Rinse solutions)
650470 Tripak Buffer Capsules (3 vials of 4, 7, 10 pH – 10 capsules per vial)
PH4-P pH 4 Buffer Solution (2 bottles – 1 pint each)
PH7-P pH 7 Buffer Solution (2 bottles – 1 pint each)
PH10-P pH 10 Buffer Solution (2 bottles – 1 pint each)
EC-84-P 84 μ S Conductivity Standard (2 bottles – 1 pint each)
EC-1413-P 1413 μ S Conductivity Standard (2 bottles – 1 pint each)
EC-12880-P 12880 μ S Conductivity Standard (2 bottles – 1 pint each)